



SAFE BABIES

INFANT EMERGENCY EVACUATION APRON



Safe Babies® Family

- Each apron has two pouches in the front and two in the back which allow for a total of four infants to be carried to safety.
- Only those staff who are physically capable of carrying the weight of four infants should be assigned to wearing the Safe Babies® apron during evacuation procedures.





Safe Babies® Family

- The Safe Babies® apron has a universal front and back orientation to help reduce confusion while putting it on during an emergency situation.





Safe Babies® Family

- Infants should be individually placed in the Safe Babies® pouches while still wrapped in their warm blankets or clothing.





Safe Babies® Family

- **WARNING** – Do not place active older infants/babies in the pouches who could possibly squirm out on their own during emergency transport.





Safe Babies® Family

- Staff members' arms are free to open doors, carry other medical supplies, or to help other staff with additional infants.





Safe Babies® Family


- The Velcro® side fasteners are to be used for quick, strong and adjustable closures.





Safe Babies® Family

- The apron is made of heavy duty, durable and washable cotton fabric.
- Machine wash in cold water using no bleach and hang dry.



**DO NOT PUT
IN DRYER!
HANG DRY
REQUIRED!**

